

Mick® Gold Fiduciary Markers



Mick® Gold Fiduciary Markers

Designed for superior image guidance during EBRT, IMRT and HDR brachytherapy procedures.

Fiducial markers in preloaded, pre-waxed needles to support visualization during radiation therapy.

Image Guidance for Radiation Therapy

Gold markers are used to optimize target visualization during external beam radiation therapy and brachytherapy procedures. Their placement in soft tissue or organs prior to therapeutic treatments help to provide for clearer identification of anatomic regions around a target treatment site. As a result, their use may aid in better dosimetric coverage of the target site.

Markers in Different Lengths and Diameters

According to specific needs and proposed treatment site, customers may choose from a variety of marker sizes. Mick® Gold Fiduciary Markers are available in 3.0 mm and 5.0 mm lengths as well as in 1.2 mm and 1.0 mm diameter. Accordingly, they come preloaded in 17 or 18G needles. All markers are cast from solid gold and provide excellent visualization under fluoroscopy, X-ray, CT, MRI, ultrasound and other imaging modalities.

Individual Packaging

Mick® Gold Fiduciary Markers are available in different packaging configurations. According to the planned treatment site and required marker quantity, they may be ordered individually packed or in 3/pouch.

Convenience through Preloaded Needles


Mick® Gold Fiduciary Markers are preloaded in pre-waxed needles, shipped sterile and are ready for immediate use thus allowing for more efficient procedures.

Mick® Gold Fiduciary Markers Characteristics

- Preloaded in a pre-waxed tipped needle
- Different diameters available: 1.0 mm or 1.2 mm
- Different lengths available: 3.0 mm or 5.0 mm

Mick® Gold Fiduciary Markers – Sizes and Packaging Units

Article No.	Needle	Marker Size	Packaging	Quantity
1421-0303	17G x 20 cm	1.2 x 3.0 mm	1/pouch	10 pouches = 10 needles
1421-0304	17G x 20 cm	1.2 x 3.0 mm	3/pouch	10 pouches = 30 needles
1421-0305	17G x 20 cm	1.2 x 5.0 mm	1/pouch	10 pouches = 10 needles
1421-0306	17G x 20 cm	1.2 x 5.0 mm	3/pouch	10 pouches = 30 needles
1421-0307	18G x 20 cm	1.0 x 3.0 mm	1/pouch	10 pouches = 10 needles
1421-0308	18G x 20 cm	1.0 x 3.0 mm	3/pouch	10 pouches = 30 needles
1421-0309	18G x 20 cm	1.0 x 5.0 mm	1/pouch	10 pouches = 10 needles
1421-0310	18G x 20 cm	1.0 x 5.0 mm	3/pouch	10 pouches = 30 needles



Mick® Gold Fiduciary Markers are each individually preloaded in a pre-waxed needle



Eckert & Ziegler

Mick Radio-Nuclear Instruments, Inc.

An Eckert & Ziegler BEBIG Company

Mick is a registered trademark of Eckert & Ziegler BEBIG s.a. and its subsidiaries.

The mentioned products are not available in all markets. Please contact your local Eckert & Ziegler BEBIG representative for more information.

Corporate Head Office:

**Eckert & Ziegler
BEBIG s.a.**
Rue Jules Bordet
Zone Industrielle C
7180 Seneffe
Belgium

Phone +32 64 520 811
Fax +32 64 520 801
info@bebig.com

Manufacturer:

CP Medical, Inc.
803 NE 25th Avenue
Portland, OR 97232
USA

Regional Sales, Marketing and Service:

Europe, Middle East, Africa,
Latin America, Asia Pacific

**Eckert & Ziegler
BEBIG s.a.**
Rue Jules Bordet
Zone Industrielle C
7180 Seneffe
Belgium

Phone +32 64 520 811
Fax +32 64 520 801
info@bebig.com

North America

**Mick Radio-Nuclear Instruments, Inc.
An Eckert & Ziegler BEBIG Company**
521 Homestead Avenue

Mount Vernon, NY 10550
USA

Phone +1 914 667 3999
Fax +1 914 665 8834
sales@micknuclear.com

www.bebig.com
www.micknuclear.com